



*Communications
for today... and tomorrow*

Tek-CARE® NC110, NC110A Nurse Call System

Submittal Package

as of 08/18/2014

www.tektone.com

Phone: (828) 524-9967
Toll-Free: (800) 327-8466
Sales: Option 2
Tech Support: Option 3

277 Industrial Park Road
Franklin, NC 28734
tektone@tektone.net
Fax: (828) 524-9968

Copyright © 2014 TekTone® Sound & Signal Mfg., Inc., All rights reserved.

No part of these publications may be copied without the express written permission of TekTone® Sound & Signal Mfg., Inc. The contents of these specification sheets are furnished for informational use only, are subject to change without notice, and should not be construed as a commitment by TekTone® Sound & Signal Mfg., Inc. TekTone® Sound & Signal Mfg., Inc. assumes no responsibility or liability for any errors or inaccuracies that may appear in this documentation.

TekTone, the TekTone logo, Tek-Alert, Tek-Bridge, Tek-Call, Tek-Care, TekCD, Tek-Check-In, Tek-Com, Tek-Entry, Tek-Guard, Tek-Micro, Tek-MMARS, Tek-Paging, Tek-Safe and Tek-Status are either registered trademarks or trademarks of TekTone® Sound & Signal Mfg., Inc. in the United States and/or other countries. All other trademarks are the property of their respective owners.



NC110-series Master Panels

IL365 Specification Sheet

Tek-CARE®NC110 Nurse Call System • Rev. 13 - 08/2013

Features

- UL® 1069 Listed
- Electronic DC buzzer
- Tone-Off and Tone-Test functions
- Emergency and Routine call signaling
- Modular design
- Protected markings

Description

TekTone®'s Tek-CARE®NC110 series Audible-Visual Nurse Call System provides audible and visual call indication without the use of intercom communications. The basic system consists of a master annunciator/control panel, room call stations, emergency call stations, dome lights, zone lights, control unit, and power supply equipment.

Emergency signals are identified at the master station by a rapid flashing of the associated station indicator light, call light, and emergency light; in addition to a rapid pulsating tone. Routine calls are indicated by the steady illumination of the station selector light, call light, and a slow pulsating tone. The tone may be canceled for normal calls by depressing the tone-off button. Emergency calls take precedence over Routine calls. Call signals can be canceled at point of origin only.

Tone-Off and Tone-Test functions are push button selected from the control panel. Indicator lights signal the presence of Emergency or Routine calls.

Specifications

Dimensions: Refer to OH200-series housings and OF200-series frames in the IL312 Housings Reference Chart

Power Requirement: 24 VAC

Connections: Screw terminal

Construction: Extruded Aluminum Panels



NC110/32 Master Panel with OF201 Frame

Listings



UL® 1069 listed; UL® File Number E71599.

Required Components

AM110-series Annunciator Add-On Panels
NC110-series Master Annunciator Panel
OH200-series Flush Housing
OF200-series Frame
SF-series Call Stations
IH151NK Power and Control Kit

Optional Components

PK124K Battery Backup Kit
PK153 Third Priority Control Unit

www.tektone.com

277 Industrial Park Road • Franklin, NC 28734 • tektone@tektone.net

Phone: (828) 524-9967 • Fax: (828) 524-9968 • Sales: Choose option 2 • Tech Support: Choose option 3
TekTone®'s quality system is registered by UL® to the ISO 9001 standard. (Reference #10001510.)



NC110A Annunciator/Control Panel

IL850 Specification Sheet

Tek-CARE®NC110 Nurse Call System • Rev. 4 - 08/2013

Features

- UL® 1069 Listed
- Electronic tone circuit with volume control
- Tone on/off and tone test functions
- Emergency and Routine call signaling
- 50-point capacity with labels for custom marking
- Spill-resistant membrane switch

Description

TekTone®'s NC110A Audible-Visual Nurse Call Annunciator/Control Panel provides for audible and visual call indication without the use of intercom communications. A basic nurse call system may consist of an NC110A annunciator/control panel, room call stations, emergency call stations, dome lights, zone lights, control unit and power supply equipment.

Emergency signals are identified at the annunciator/control panel by a rapid flashing of the associated station indicator LED, call LED and emergency LED; in addition to a rapid pulsating tone. Routine calls are indicated by the steady illumination of the station selector LED, call LED, and a slow pulsating tone. The tone may be canceled for routine calls by depressing the tone on/off button. Emergency calls take precedence over Routine calls. Emergency call signals can be canceled at point of origin only.

Tone on/off and tone test functions are button selectable from the control panel, and three levels of tone volume are button selectable with LED indication. Indicator LEDs signal the presence of Emergency or Routine calls. The NC110A has a capacity of 50 selector points. All selector LEDs, indicator LEDs and control buttons are included on a membrane-type switch mounted on a telephone-sized plastic housing. A normally open isolated relay contact output for ancillary off-board call handling is also included.

Listings



UL® 1069 listed; UL® File Number E71599.



NC110A Annunciator/Control Panel

Specifications

Dimensions:	8.5" × 8.5" × 3" (216 mm × 216 mm × 76 mm)
Power Requirement:	24 VDC, 650 mA
Connections:	Keyed plug-in
Construction:	Flame-retardant (UL®94V-0) plastic with membrane switch
Mounting:	3-gang ring for connector plate

Required Components

SF-series Call Stations
IH151NK Power & Control Kit

Optional Components

PK124K Battery Backup Kit
PK153 Third Priority Control Unit

www.tektone.com

277 Industrial Park Road • Franklin, NC 28734 • tektone@tektone.net

Phone: (828) 524-9967 • Fax: (828) 524-9968 • Sales: Choose option 2 • Tech Support: Choose option 3
TekTone®'s quality system is registered by UL® to the ISO 9001 standard. (Reference #10001510.)



IH300-series Desk Cabinets

IL502 Specification Sheet

CM800, DS100, and Tek-CARE®NC110 & NC150 Systems • Rev. 11 - 08/2013



IH302 Housing



IH303 Housing

Description

IH300-series housings are desk cabinets for 2 to 6 TekTone® 4" wide aluminum panels. See chart for reference. IH300-series cabinets are constructed of black laminate covered wood.

Required Components

Housings are designed for use with Tek-CARE® NC110 and Tek-CARE® NC150 nurse call panels, CM800 emergency call panels and DS100 door sentry & emergency alarm panels.

Specifications

Model #	Number of Panels	Height	Width	Depth
IH302	2	6"	9.5"	17"
IH303	3	6"	13.5"	17"
IH304	4	6"	17.5"	17"
IH305	5	6"	21.5"	17"
IH306	6	6"	25.5"	17"

www.tektone.com

277 Industrial Park Road • Franklin, NC 28734 • tektone@tektone.net

Phone: (828) 524-9967 • Fax: (828) 524-9968 • Sales: Choose option 2 • Tech Support: Choose option 3
TekTone®'s quality system is registered by UL® to the ISO 9001 standard. (Reference #10001510.)



Housings Reference Chart

IL312
Rev. 20 - 06/2013

NOTES:

1. For SF-series Emergency Switches, use single-gang ring, or single-gang ring and double-gang box IH357/IH358.
2. For LI-series Dome Lights, use double-gang ring (or double-gang ring and box).
3. For AM- and CM-series, see specification requirements.
4. For IR300-series Nurse Call Stations, use Steel City #3GC and #H3BD or exact equivalent.
5. For IR150B-series, use Steel City #3GC and #H3BD. If space is limited, use ARC-CO #2GBD and #3GC.
6. For IR150SS Station, use OH600 flush housing mounted horizontally.
7. The last two digits of frame and housing model numbers indicate the number of gangs (i.e., panel capacity). Exceptions: IH101, IH103, OH600, OH601, IH400A.
8. Attach the IH103 mounting plate to a finished wall, single-gang box, double-gang box or square masonry box.

Model #	IH101 backbox (flush)	IH103 plate (flush)	IH400A housing (desk)
IR100B	x	x	x
IR102C	x	x	x
IR103E	x	x	x
IR104E	x	x	x
IR105E	x	x	x
IR203E		x	x
IR204E		x	x
IR205E		x	x
IR207B	use standard 2-gang electrical box		
IR704BD	x	x	x
OR500 series	x	x	
OR600 series	use standard 1-gang electrical box		
SM201B SM201K/P	use standard 1-gang electrical box		x

Dimension and Finish Chart—Housings, Frames, Desk Cabinets

Flush Mount Housings					
Model#	Overall Width	Wall Cutout Width	Overall Height	Wall Cutout Height	Overall Depth
IH101	4"	4-1/4"	5-1/2"	5-3/4"	1-3/4"
IH103	4-9/16"	n/a	5-7/8"	n/a	1/16"
OH150	5-1/4"	4-1/2"	12-5/8"	11-5/8"	3-3/4"
OH191	4-3/8"	4-1/2"	17-1/8"	17-1/4"	3-1/2"
OH192	8-3/8"	8-1/2"	17-1/8"	17-1/4"	3-1/2"
OH193	12-3/8"	12-1/2"	17-1/8"	17-1/4"	3-1/2"
OH194	16-3/8"	16-1/2"	17-1/8"	17-1/4"	3-1/2"
OH195	20-3/8"	20-1/2"	17-1/8"	17-1/4"	3-1/2"
OH201	5-3/8"	4-1/2"	16-7/8"	16-1/4"	3-7/8"
OH202	9-3/8"	8-1/2"	16-7/8"	16-1/4"	3-7/8"
OH203	13-3/8"	12-1/2"	16-7/8"	16-1/4"	3-7/8"
OH600	4"	4-1/4"	10-3/16"	10-3/8"	1-13/16"
OH601	4"	4-1/4"	16-1/16"	16-1/4"	1-13/16"

Flush Mount Frames		
Model #	Width	Wall Cutout
OF191	6-1/4"	-
OF192	10-1/4"	-
OF193	14-1/4"	-
OF194	18-1/4"	-
OF195	22-1/4"	-
Height: 19" (overall)		
Depth: 7/8"		
OF201	6-1/4"	4-1/2"
OF202	10-1/4"	8-1/2"
OF203	14-1/4"	12-1/2"
Height: 17-3/4" (overall)		
Cutout: 15-7/8"		
Depth: 2-1/2"		

Desk Cabinets		
Model #	Width	Depth
*IH302	9-1/2"	6"
*IH303	13-1/2"	6"
*IH304	17-1/2"	6"
IH400A (white)	5"	7"
* Cabinet Assembly Only. IH300 series Height = 17"		

Surface Mount Frames			
Model #	Width	Height	Depth
OF190S	4-1/4"	19-3/4"	1-1/2"
OF191S	6-1/4"	19"	2-1/4"
OF192S	10-1/4"	19"	2-1/4"
OF193S	14-1/4"	19"	2-1/4"
OF194S	18-1/4"	19"	2-1/4"
OF195S	22-1/4"	19"	2-1/4"

Surface Mount Housings	
Model #	Width
OH301	6-1/8"
OH302	10-1/8"
Height: 17-3/4"	
Depth: 3-5/8"	

www.tektone.com

277 Industrial Park Road • Franklin, NC 28734 • tektone@tektone.net

Phone: (828) 524-9967 • Fax: (828) 524-9968 • Sales: Choose option 2 • Tech Support: Choose option 3
TekTone®'s quality system is registered by UL® to the ISO 9001 standard. (Reference #10001510.)



SF100C Single Bed Room Station

IL362 Specification Sheet

Tek-CARE® NC110 Nurse Call System • Rev. 14 - 08/2013

Features

- UL® 1069 Listed
- Simple, push button calling using plug-in call cord
- Call-placed indicator and cancel button
- Call automatically registers if call cord is accidentally removed
- Corridor zone light capability

Description

The SF100C Room Station includes a jack allowing the attachment of an SF301-series call cord. A call is registered at the nurse station, corridor light and corridor zone light (if installed) when the call cord button is used, or if the call cord is accidentally removed from the jack on the panel. The station includes a call-placed indicator that lights when a call is placed and a cancel button mounted on an attractive, flame-retardant (UL® 94V-0) plastic faceplate.

Listings



UL® 1069 listed; UL® File Number E71599.



SF100C Single Bed Room Station

Specifications

Dimensions:	Height:	4.5"	(114 mm)
	Width:	2.75"	(69 mm)
	Depth:	1.75"	(44 mm)
	Projects:	0.4375"	(11 mm)
Construction:	Attractive, flame retardant (UL® 94V-0) plastic faceplate. Switches and electronic components are circuit board mounted.		
Connections:	Plug-in connectors.		
Housing:	Single gang ring, or single gang ring with double gang box.		

Required Components

Tek-CARE® NC110 Nurse Call System
SF301-series Call Cord

www.tektone.com

277 Industrial Park Road • Franklin, NC 28734 • tektone@tektone.net

Phone: (828) 524-9967 • Fax: (828) 524-9968 • Sales: Choose option 2 • Tech Support: Choose option 3

TekTone®'s quality system is registered by UL® to the ISO 9001 standard. (Reference #10001510.)



SF101C Audible-Visual Room Station

IL366 Specification Sheet

Tek-CARE® NC110 Nurse Call System • Rev. 12 - 08/2013

Features

- UL® 1069 Listed
- Simple, push-button calling
- Call-placed indicator
- Cancel push button
- Call activates lamp on annunciator panel and corridor light
- Flame retardant (UL® 94V-0) plastic faceplate
- Corridor zone light capability

Description

The SF101C Room Station incorporates push-button operation of call and cancel functions. A call is registered at the annunciator panel and on a corridor light or corridor zone light (if installed) when the call switch with integral call-placed LED is pushed. The call-placed LED will also illuminate. Calls are canceled by simply depressing the cancel button. The call-placed LED will also be extinguished. All station components are mounted on a rugged UL® 94V-0 plastic faceplate.

Listings



UL® 1069 listed; UL® File Number E71599.



SF101C Audible-Visual Room Station

Specifications

Dimensions:	Height:	4.50"	(114 mm)
	Width:	2.75"	(69 mm)
	Depth:	1.75"	(44 mm)
	Projects:	0.4375"	(11 mm)
		from mounting surface	
Construction:	Flame retardant (UL® 94V-0) plastic faceplate. Switches and electronic components are circuit board mounted.		
Connections:	Plug-in connectors.		
Housing:	Single gang ring, or single gang ring with double gang box.		

Required Components

Tek-CARE® NC110 Nurse Call System

www.tektone.com

277 Industrial Park Road • Franklin, NC 28734 • tektone@tektone.net

Phone: (828) 524-9967 • Fax: (828) 524-9968 • Sales: Choose option 2 • Tech Support: Choose option 3

TekTone®'s quality system is registered by UL® to the ISO 9001 standard. (Reference #10001510.)



SF102 Dual Bed Room Station

IL639 Specification Sheet

Tek-CARE® NC110 Nurse Call System • Rev. 6 - 08/2013

Features

- UL® 1069 Listed
- Simple, push button calling using plug-in call cords
- Combined call-placed indicator and cancel button
- Call automatically registers if either call cord is accidentally removed
- Attractive, flame retardant ABS plastic faceplate
- Corridor zone light capability

Description

The SF102 Room Station includes two jacks allowing the attachment of two SF301-series call cords. A call is registered at the nurse station, corridor light and corridor zone light (if installed) when either cord button is depressed or if either call cord is accidentally removed from the jack on the SF102.

The station includes a lighted call-placed switch that lights when a call is placed. The call-placed switch also provides the system cancel function. All of the electronics are mounted on an attractive, flame retardant ABS plastic faceplate.

Listings



UL® 1069 listed; UL® File Number E71599.



SF102 Dual Bed Room Station

Specifications

Dimensions:	Height:	4.5"	(114 mm)
	Width:	2.75"	(69 mm)
	Depth:	1.75"	(44 mm)
	Projects:	0.4375"	(11 mm)
Construction:	Attractive, flame retardant ABS plastic faceplate. Switch and electronic components are circuit board mounted.		
Connections:	Plug-in connectors.		
Housing:	Single gang ring, or single gang ring with double gang box.		

Required Components

Tek-CARE® NC110 Nurse Call System
SF301-series Call Cords

www.tektone.com

277 Industrial Park Road • Franklin, NC 28734 • tektone@tektone.net

Phone: (828) 524-9967 • Fax: (828) 524-9968 • Sales: Choose option 2 • Tech Support: Choose option 3
TekTone®'s quality system is registered by UL® to the ISO 9001 standard. (Reference #10001510.)



SF301A, SF302 Call Cords IL338 Specification Sheet

Tek-CARE® Nurse Call Systems • Rev. 18 - 06/2013

Features

- UL® 1069 Listed
- Oxygen safe
- Impact resistant and waterproof
- Small, lightweight pendant fits comfortably in hand
- Bed clips included

Description

The TekTone® SF301A (single) or SF302 (double) call cords are for use with the Tek-CARE® NC110N, NC150N, NC200N and NC300™II Nurse Call Systems, and the Tek-CARE®500 Wireless Emergency Call System, where ¼" phono jacks are used. Ten-foot cords are also available as SF301A/10 and SF302/10. The SF301A and SF302 are oxygen safe call cords certified safe for use by all patients receiving oxygen therapy. These waterproof and impact resistant call cords have completely sealed push buttons.



SF301A, SF302 Call Cords

Listings



UL® 1069 Listed; UL® File Number E71599.

Specifications

Connection: ¼" (6.35 mm) plug-in phono jack
 Construction: SPST-1A switch
 Integral brass, fine silver contacts
 1 AMP, 115 watt maximum, AC
 Dimensions: 7' (2.13 m) long polyvinyl covered cord

Required Components

Tek-CARE®500 Wireless Emergency Call System
with SF530 Room Station

or

Tek-CARE® NC300™II Nurse Call System
with SF380A Patient Stations

or

Tek-CARE® NC200N or NC150N Nurse Call Systems
with IR151B or IR152B Patient Stations

or

Tek-CARE® NC110N Nurse Call System
with SF100C or SF102 Patient Stations

www.tektone.com

277 Industrial Park Road • Franklin, NC 28734 • tektone@tektone.net

Phone: (828) 524-9967 • Fax: (828) 524-9968 • Sales: Choose option 2 • Tech Support: Choose option 3

TekTone®'s quality system is registered by UL® to the ISO 9001 standard. (Reference #10001510.)



SF301G, SF302G Call Cords IL434 Specification Sheet

Tek-CARE® Nurse Call Systems • Rev. 15 - 06/2013

Features

- UL® 1069 Listed
- Sensitive touch activation
- Oxygen tent safe
- Ideal for geriatric use
- Momentary, air-activated switch
- Adjustable bed clip(s)

Description

TekTone® SF301G and SF302G Call Cords are for use with the Tek-CARE® NC110N, NC150N, NC200N and NC300II Nurse Call Systems, where ¼" phono jacks are used. Patient stations include SF100C, IR151B, IR152B, SF380C, SF380P and SF530.

SF301G and SF302G Call Cords are specifically designed for use in geriatric and oxygen environments, but may also be used in general applications. To place a call, the patient applies slight pressure to the squeeze bulb using any part of the body. The force of air through the vinyl tubing closes a switch in the plug, thus placing a call. An adjustable bed clip is supplied for the SF301G and is attached to the call cord. Two adjustable bed clips are supplied with the SF302G.

Listings



UL® 1069 Listed; UL® File Number E71599.

Specifications

Switch:	Momentary contact, air activated		
Plug:	Stainless steel		
Cord:	Vinyl tubing		
Electrical:	Rated current 1/3 amp, 24 VAC-VDC		
Dimensions:	Cord:	6' long	(1.83 m)
	Plug:	1.25" long	(32 mm)
	Bulb:	1.5" diameter	(38 mm)



SF301G Call Cord

Required Components

Tek-CARE®500 Wireless Emergency Call System
with SF530 Room Station

or

Tek-CARE® NC300™II Nurse Call System
with SF380A Patient Stations

or

Tek-CARE® NC200N or NC150N Nurse Call Systems
with IR151B or IR152B Patient Stations

or

Tek-CARE® NC110N Nurse Call System
with SF100C or SF102 Patient Stations

www.tektone.com

277 Industrial Park Road • Franklin, NC 28734 • tektone@tektone.net

Phone: (828) 524-9967 • Fax: (828) 524-9968 • Sales: Choose option 2 • Tech Support: Choose option 3

TekTone®'s quality system is registered by UL® to the ISO 9001 standard. (Reference #10001510.)



SF155B Emergency Station

IL513 Specification Sheet

Tek-CARE® NC110, NC150, NC200 Nurse Call Systems • Rev. 11 - 08/2013

Features

- UL® 1069 Listed
- Simple and flexible operation
- Call-placed indicator
- Call cancel at point of origin
- Flame-retardant (UL® 94V-0) plastic faceplate
- 1.1" wide pull down lever

Description

The SF155B Emergency Station is for use in any location where activation of an emergency signal is required. All electronic components are mounted on a circuit board. The specially designed and mounted LED is completely sealed in plastic.

Due to the flexibility of its design, it may be used as a pull-cord emergency station, or a pull-down activated emergency station. The pull cord supplied may be used for emergency call activation, or if desired, the large 1.1" red pull-down lever may be used.

Listings



UL® 1069 listed; UL® File Number E71599.



SF155B Emergency Station

Specifications

Dimensions:	Height:	4.5"	(114 mm)
	Width:	2.8125"	(72 mm)
	Depth:	1.625"	(41 mm)
	Projects:	0.6875"	(17 mm)
		from mounting surface (including pull-down lever)	
Construction:	Flame-retardant (UL® 94V-0) plastic faceplate. Switch and electronic components are circuit board mounted.		
Connections:	Plug-in connectors.		
Housing:	Single gang ring, or single gang ring with double gang box.		

Required Components

Tek-CARE® NC110, NC150 or NC200 Nurse Call System

www.tektone.com

277 Industrial Park Road • Franklin, NC 28734 • tektone@tektone.net

Phone: (828) 524-9967 • Fax: (828) 524-9968 • Sales: Choose option 2 • Tech Support: Choose option 3

TekTone®'s quality system is registered by UL® to the ISO 9001 standard. (Reference #10001510.)



SF156B Code Station

IL648 Specification Sheet

Tek-CARE® NC110, NC150, NC200 Nurse Call Systems • Rev. 8 - 08/2013

Features

- UL® 1069 Listed
- Simple, reliable operation
- Call-placed indicator
- Call activates light on annunciator panel, corridor light and/or zone light.
- Corridor zone light capability
- Mounted on attractive flame retardant (UL® 94V-0) plastic panel

Description

The SF156B Code Station is used to operate all system code signals (except supervised code blue applications). It may be used independently or connected with other call stations.

Placing a code call is accomplished by pulling the slide-switch located at the center of the panel to the down position. Returning the slide-switch to the up position cancels the call.

Listings



UL® 1069 listed; UL® File Number E71599.



SF156B Code Station

Specifications

Dimensions:	Height:	4.5"	(114 mm)
	Width:	2.75"	(70 mm)
	Depth:	0.75"	(20 mm)
Construction:	Reliable components mounted on an attractive flame retardant (UL® 94V-0) plastic panel.		
Connections:	Plug-in wiring pigtail.		
Housing:	Single gang ring, or single gang ring with double gang box.		

Required Components

Tek-CARE® NC110, NC150 or NC200 Nurse Call System
PK153 Third Priority Control Unit

www.tektone.com

277 Industrial Park Road • Franklin, NC 28734 • tektone@tektone.net

Phone: (828) 524-9967 • Fax: (828) 524-9968 • Sales: Choose option 2 • Tech Support: Choose option 3
TekTone®'s quality system is registered by UL® to the ISO 9001 standard. (Reference #10001510.)



SF337C Emergency Shower Switch

IL629 Specification Sheet

Tek-CARE® Nurse Call Systems • Rev. 12 - 11/2012

Features

- UL® 1069 Listed
- Waterproof if installed per manufacturer's instructions
- Call-placed indicator
- Call cancel at point of origin
- Attractive, flame-retardant (UL® 94V-0) plastic faceplate
- 1.1" (28mm) wide pull-down lever
- Washable TPU plastic pull cord kit is included

Description

The SF337C Emergency Shower Switch is for any location where activation of an emergency signal is required. The design is such that it may be located anywhere in a shower, bath or tub room, regardless of environmental conditions. The combination of rubber-grommetted mounting screws, rubber mounting gasket, and completely sealed front faceplate make it waterproof (if properly installed). Electronic components are mounted on a circuit board, including an industry-proven reliable reed switch for circuit initiation. The activating magnet is mounted on the large pull-down lever, which is completely sealed in plastic. The indicator LED is also completely sealed.

Due to the flexibility of its design, the SF337C may be used as a pull cord or as a pull down activated emergency switch. The pull cord supplied may be used for emergency call activation, or if desired, the large 1.1" (28 mm) red pull-down lever may be used. Dome light and system annunciation vary according to the manner in which the SF337C is connected.

Listings



UL® 1069 listed; UL® File Number E71599.



SF337C Emergency Shower Switch
(cord not shown)

Specifications

Dimensions:	Height:	4.5"	(114 mm)
	Width:	2.8125"	(72 mm)
	Depth:	1.625"	(41 mm)
	Projects:	0.6875"	(17 mm)
			(including pull down lever)
Construction:	Flame-retardant (UL® 94V-0) plastic		
Connections:	Plug-in type.		
Housing:	Single gang ring with square box.		

Applicable Systems

- Tek-CARE®400 P5 Nurse Call System
- Tek-CARE®NC300II Nurse Call System
- Tek-CARE® NC200, NC150, NC110 Nurse Call Systems
- Tek-CARE®500 Wireless Emergency Call System

www.tektone.com

277 Industrial Park Road • Franklin, NC 28734 • tektone@tektone.net

Phone: (828) 524-9967 • Fax: (828) 524-9968 • Sales: Choose option 2 • Tech Support: Choose option 3
TekTone®'s quality system is registered by UL® to the ISO 9001 standard. (Reference #10001510.)



SF339 Emergency Switch

IL744 Specification Sheet

Tek-CARE® Nurse Call Systems • Rev. 5 - 11/2012

Features

- UL® 1069 Listed
- Simple and flexible operation
- Call-placed indicator
- Call cancel at point of origin
- Flame-retardant (UL® 94V-0) plastic faceplate
- 0.9" diameter (23 mm) red mushroom-style button

Description

The SF339 Emergency Switch includes a locking mushroom-style switch, and is for use in any location where activation of an emergency signal is required. All electronic components are mounted on a circuit board, including a reliable industry-proven DPDT switch for circuit initiation. The DPDT switch and the indicator LED are completely sealed in plastic.

The SF339 may be used as an emergency switch or as a staff-assist emergency call switch. The large mushroom-style button is easily activated, so that calls may be placed using an elbow, shoulder, etc., when using the hands is inconvenient or inappropriate. Calls are cancelled by turning the mushroom-style button clockwise.

Dome light indication and system annunciation vary according to the manner in which the SF339 is installed.

Specifications

Dimensions:	Height:	4.5"	(114.3 mm)
	Width:	2.8125"	(71.4 mm)
	Depth:	1"	(25.4 mm)
	Projects:	1"	(25.4 mm)
Construction:	Flame-retardant (UL® 94V-0) plastic faceplate. DPDT switch and electronic components are circuit board mounted.		
Connections:	Plug-in connectors.		
Housing:	Single gang ring (IH357) with square back box (IH358) or exact equal.		



SF339 Emergency Switch

Listings



UL® 1069 listed; UL® File Number E71599.

Applicable Systems

- Tek-CARE®400 P5 Nurse Call System
- Tek-CARE®NC300II Nurse Call System
- Tek-CARE® NC200, NC150, NC110 Nurse Call Systems
- Tek-CARE®500 Wireless Emergency Call System

www.tektone.com

277 Industrial Park Road • Franklin, NC 28734 • tektone@tektone.net

Phone: (828) 524-9967 • Fax: (828) 524-9968 • Sales: Choose option 2 • Tech Support: Choose option 3

TekTone®'s quality system is registered by UL® to the ISO 9001 standard. (Reference #10001510.)



LI150B Duty Station

IL654 Specification Sheet

Tek-CARE® NC110, NC150, NC200 Nurse Call Systems • Rev. 5 - 08/2013

Features

- UL® 1069 Listed
- Simple, reliable operation
- Call-placed indicator
- Buzzer for audible annunciation
- Mounted on rugged, flame-retardant (UL® 94V-0) plastic panel

Description

The LI150B Duty Station provides an indication of incoming calls in work areas away from the nurse call master station.

Normal calls are indicated both visibly by steady illumination of the solid-state LED light, and audibly by a slow pulsating tone. Emergency calls take precedence and are indicated by flashing of the LED light and a rapid pulsating tone. A tone-off switch permits silencing of the tone signal for routine level calls only.

Listings



UL® 1069 listed; UL® File Number E71599.



LI150B Duty Station

Specifications

Dimensions:	Height:	4.5"	(114 mm)
	Width:	2.75"	(70 mm)
	Depth:	1"	(25 mm)
	Projects:	0.5"	(13 mm)
		from mounting surface	
Construction:	Reliable components mounted on a flame-retardant (UL® 94V-0) plastic faceplate		
Connections:	Screw terminals		
Housing:	Single gang ring, or single gang ring with double gang box.		

Required Components

Tek-CARE® NC110, NC150 or NC200 Nurse Call System
IR150-series or SF100-series Patient Stations

www.tektone.com

277 Industrial Park Road • Franklin, NC 28734 • tektone@tektone.net

Phone: (828) 524-9967 • Fax: (828) 524-9968 • Sales: Choose option 2 • Tech Support: Choose option 3
TekTone®'s quality system is registered by UL® to the ISO 9001 standard. (Reference #10001510.)



LI381 Corridor Light & LI382 Corridor Zone Light IL518 Specification Sheet

Tek-CARE® Nurse Call & CM800 Emergency Call Systems • Rev. 18 - 07/2013

Features

- UL® 1069 Listed
- Snap-in lens
- Easy bulb replacement
- Attractive flame-retardant HIPS30 plastic panel
- Excellent visibility

Description

The LI381 Corridor Light and LI382 Corridor Zone Light provide visual indication of calls when used with a nurse call system, emergency call system, or similar application. Emergency calls are indicated by the rapid flashing of the LI381 Corridor Light associated with the calling station. Normal calls produce a steady illumination.

The LI382 Corridor Zone Light indicates an emergency call by the rapid flashing of the red-covered bulb associated with the zone or area from which an emergency call has been placed. Normal calls are indicated by a steady illumination of the white bulb. If an emergency and normal call are placed in the same zone at the same time, the red covered bulb connected to the emergency station will flash rapidly, while the white bulb associated with the normal call will maintain a steady illumination.

The standard mounting plate accommodates both models and features an easily removable white flame-retardant translucent polystyrene plastic lens for bulb replacement.

Listings



UL® 1069 listed; UL® File Number E71599.



LI381 Corridor Light & LI382 Corridor Zone Light

Specifications

Dimensions:	Height:	4.5625"	(115 mm)
	Width:	4.5625"	(115 mm)
	Depth:	3.3"	(84 mm)
Bulbs:	#1829 28V 70 mA		
Power:	LI381:	24 VDC, 51mA	
	LI382:	24 VDC, 102mA	
Finish:	White flame-retardant translucent polystyrene plastic lens mounted on a flame-retardant (UL®94V-0) plastic panel.		
Housing:	Fits either standard dual gang box and ring or single gang box and ring.		

Required Components

Tek-CARE® NC110, NC150, NC200 or NC300II Nurse Call System

-or-

CM800 Emergency Call System

www.tektone.com

277 Industrial Park Road • Franklin, NC 28734 • tektone@tektone.net

Phone: (828) 524-9967 • Fax: (828) 524-9968 • Sales: Choose option 2 • Tech Support: Choose option 3

TekTone®'s quality system is registered by UL® to the ISO 9001 standard. (Reference #10001510.)



LI382LED Corridor Zone Light IL898 Specification Sheet

Tek-CARE® Nurse Call & CM800 Emergency Call Systems • Rev. 2 - 07/2013

Features

- UL® 1069 Listed
- Long lasting LEDs
- No bulb replacement required
- Plug-in connector supplied
- Attractive flame-retardant plastic panel

Description


The LI382LED Corridor Zone Light provides visual indication of calls when used with a nurse call system, emergency call system, or similar application.

The LI382LED Corridor Zone Light indicates an emergency call by the rapid flashing of the two red LED bulbs associated with the zone or area from which an emergency call has been placed. Normal calls are indicated by a steady illumination of the two white LEDs. If an emergency and normal call are placed in the same zone at the same time, the red LEDs connected to the emergency station will flash rapidly, while the white LEDs associated with the normal call will maintain a steady illumination.

The standard mounting plate features both one- and two-gang mounting options.

Specifications

Dimensions:	Height:	4.8"	(122 mm)
	Width:	5.0"	(127 mm)
	Depth:	2.755"	(70 mm)
LEDs:	2V 20mA per LED		
Power:	24 VDC		
Finish:	Base is flame-retardant (UL® 94HB) plastic. Lens is flame-retardant translucent polycarbonate plastic.		
Housing:	Fits standard one- or two-gang box with one- or two-gang ring.		



LI382LED Corridor Zone Light

Listings



UL® 1069 listed; UL® File Number E71599.

Required Components

Tek-CARE® NC110, NC150, NC200 or NC300II Nurse Call System

-or-

CM800 Emergency Call System

www.tektone.com

277 Industrial Park Road • Franklin, NC 28734 • tektone@tektone.net

Phone: (828) 524-9967 • Fax: (828) 524-9968 • Sales: Choose option 2 • Tech Support: Choose option 3

TekTone®'s quality system is registered by UL® to the ISO 9001 standard. (Reference #10001510.)



IH151NK Power & Control Kit IL367 Specification Sheet

Rev. 14 - 04/2013

Features

- UL® 1069 Listed power & control unit
- UL® Listed Class 2 transformer.
- Provides steady light and slow tone for normal call; flashing light and rapid tone for emergency call.
- Routine and emergency commons incorporate short-protection circuitry.

Description

The IH151NK Power & Control Kit includes a UL® 1069 Listed PK152 Power & Control Unit and an SS106 UL® Listed Class 2 Transformer, preinstalled in a 19-gauge cold-rolled steel housing. Knockouts are included on all sides except the lid, to allow connection of electrical conduit from any direction. The lid has a knockout which can be used to install an optional key lock.

The PK152 Power & Control Unit contains the electronic circuitry necessary to detect the presence of normal and emergency calls, and to generate audible and visual control signals. It provides steady or intermittent voltages to operate visual signaling devices. The PK152 supports up to 1 amp of current load from field devices, and supports interconnection of an auxiliary power supply (PK800A) to increase the output current load to 3 amps. Protective circuitry has been incorporated to protect routine and emergency call commons from shorts.



IH151NK Power & Control Kit

Specifications

Dimensions:	Height:	11.75"	(300 mm)
	Width:	11.75"	(300 mm)
	Depth:	3.75"	(96 mm)
Construction:	19 gauge, cold-rolled steel		
Transformer:	24 VAC, 30 VA (SS106)		
Connections:	Pigtails (on PK152)		
Temperature:	0°C–50°C.		
Mounting:	Surface-mounted within 50' of master station		

Specifications are the minimum requirement for system operation. Some systems may require additional equipment. Contact your factory representative for details.

Required Components

Tek-CARE® NC110, NC150 or NC200 Nurse Call System
-or- Tek-SAFE® Area of Rescue Assistance System

For higher current applications add:

- PK800A 24 VDC Auxiliary Power Supply
- FZ151 In-line 4A Fuse Assembly
- Edwards #88-100 24 VAC, 100VA Transformer, or equivalent

Optional Components

Mier BW3000 Offset Cam Lock, or equivalent

www.tektone.com

277 Industrial Park Road • Franklin, NC 28734 • tektone@tektone.net

Phone: (828) 524-9967 • Fax: (828) 524-9968 • Sales: Choose option 2 • Tech Support: Choose option 3
TekTone®'s quality system is registered by UL® to the ISO 9001 standard. (Reference #10001510.)



SS102A and SS106 Transformers PO001 Electric Door Opener IL127 Specification Sheet

Rev. 14 - 04/2013



SS102A



SS106

SS102A and SS106 Transformers



PO001 Electric Door Opener

Transformer Features

TekTone®'s SS102A and SS106 Class 2 Transformers offer compact, easy-to-install low voltage AC power sources for use with various TekTone® intercom systems. Both transformers are UL® Listed Class 2 transformers and will safely take themselves offline when overheated.

The transformers are constructed for easy surface installation or in a standard 1/2" knockout, and may be used in plastic or metal boxes. Primary wire connections are made using the pre-stripped pigtailed and ground wire; secondary connections are made using screw terminals.

Specifications

SS102A	120 VAC primary, 16 VAC 10 VA secondary
SS106	120 VAC primary, 8 VAC 20 VA 16/24 VAC 30 VA secondary

Electric Door Opener Features

TekTone®'s PO001 Electric Door Opener is for use in unlocking and locking doors from a variety of TekTone® remote stations. The door opener can be installed in place of the door lock, or on a door jamb with various door lock sets.

The PO001 is a normally locked type with heavy duty coils, and operates on 8-16 VAC or 3-6 VDC. The unit features a flat face with satin zinc finish and solid brass locking latch for long wear. The PO001 is reversible for mortise-type locks and latches up to 1/2" width and throw, and is ideal for most narrow jamb mortise-type installations. Depending on the specific job requirements, the PO001 can be reversed for right or left-hand door installation.

Dimensions

Faceplate:	1.25" × 5.875"	(31.8 mm × 149 mm)
Case width:	0.9375"	(24 mm)
Case depth:	1.4375"	(37 mm)
Latch opening:	1.375"	(35 mm)

www.tektone.com

277 Industrial Park Road • Franklin, NC 28734 • tektone@tektone.net

Phone: (828) 524-9967 • Fax: (828) 524-9968 • Sales: Choose option 2 • Tech Support: Choose option 3
TekTone®'s quality system is registered by UL® to the ISO 9001 standard. (Reference #10001510.)



PK153 Third Priority Control Unit

IL473 Specification Sheet

Tek-CARE® NC110, NC150, NC200 Nurse Call Systems • Rev. 11 - 08/2013

Features

- UL® 1069 Listed
- High-quality integrated circuitry and electronic components
- Rugged extruded aluminum case
- Provides flashing light and rapid tone for emergency calls

Description

The PK153 Third Priority Control Unit provides the control circuits to operate a higher priority call signal, such as a smoke detector alarm or a code call, for TekTone®'s Tek-CARE® NC110, NC150, NC200 nurse call systems. When used in conjunction with the IH151NK Power & Control Kit, the unit provides the additional circuitry for a third level of call priority.

Listings



UL® 1069 listed; UL® File Number E71599.

Specifications

Dimensions:	Height:	9.5"	(240 mm)
	Width:	2.375"	(60 mm)
	Depth:	2.25"	(57 mm)
Construction:	Electronic components mounted on printed circuit boards and housed in a rugged extruded aluminum case		
Connections:	Screw terminals		
Temperature Range:	0°C–50°C		
Mounting:	Surface mounted (within 50' of Nurse Call Master)		
Power Requirement:	24 VDC (provided by IH151NK)		



PK153 Third Priority Control Unit

Required Components

Tek-CARE® NC110, NC150 or NC200 Nurse Call System
IH151NK Power and Control Kit

Note: Specifications are the minimum requirement for system operation. Some systems may require additional equipment. Contact your factory representative.

www.tektone.com

277 Industrial Park Road • Franklin, NC 28734 • tektone@tektone.net

Phone: (828) 524-9967 • Fax: (828) 524-9968 • Sales: Choose option 2 • Tech Support: Choose option 3
TekTone®'s quality system is registered by UL® to the ISO 9001 standard. (Reference #10001510.)



PK800A Auxiliary Power Supply

IL664 Specification Sheet

Tek-CARE® NC110, NC150, NC200 Nurse Call Systems • Rev. 9 - 07/2014

Features

- UL® 1069 Listed
- Compact design
- Protected from shorts at the outputs

Description

TekTone®'s PK800A Auxiliary Power Supply provides additional DC power for the Tek-CARE® NC110, NC150 and NC200 Nurse Call Systems, allowing for extended light capacity. Both of the outputs, +24VDC and +5VDC, are protected from shorts that might occur at the output. When the shorts are removed, the power supply recovers and operates normally. A 4 amp fuse protects the input of the power supply.

The PK800A Auxiliary Power Supply maintains its output for approximately 2 seconds after removal of its input power. This feature prevents the loss of system functionality during periods of brief power interruptions.



PK800A Auxiliary Power Supply

Listings



UL® 1069 listed; UL® File Number E71599.

Specifications

Dimensions:	Height:	10.625"	(270 mm)
	Width:	3.5"	(88.9 mm)
	Depth:	2.5"	(64 mm)
Construction:	Extruded Aluminum case. Internal circuit board.		
Connections:	Screw Terminals		
Fuse:	4 amp		
Power:	24VAC / 100VA / 60Hz		

Required Components

Tek-CARE® NC110, NC150 or NC200 Nurse Call System
 IH151NK Power and Control Kit
 24VAC, 100VA transformer supplied by others

www.tektone.com

277 Industrial Park Road • Franklin, NC 28734 • tektone@tektone.net

Phone: (828) 524-9967 • Fax: (828) 524-9968 • Sales: Choose option 2 • Tech Support: Choose option 3

TekTone®'s quality system is registered by UL® to the ISO 9001 standard. (Reference #10001510.)



PK124K Battery Backup Kit

IL818 Specification Sheet

Rev. 5 - 04/2013

Features

- Low voltage operation
- Reliable alternate power supply
- Integrated circuit design
- Easy installation

Description

The PK124K Battery Backup Kit provides continuous battery charging so that if power to system equipment is interrupted, standby operating power is provided for up to approximately 6 hours. (Time may vary depending on system size and usage.)

The PK124K kit includes a transformer, a battery charger and two batteries—all preinstalled in a 19-gauge cold-rolled steel housing. Knockouts are included on all sides except the lid, to allow connection of electrical conduit from any direction. In addition, the lid has a knockout which can be used to install an optional key lock. The enclosed single-gang electrical junction box is galvanized steel with a stainless steel cover plate.

TekTone® systems supported by the PK124K include the DS100 Door Sentry & Emergency Alarm System, the CM800 Emergency Call System, Tek-ENTRY® Apartment Entry Systems, and the Tek-SAFE® Area of Rescue Assistance System.

Specifications

Dimensions:	Height:	11.75"	(300 mm)
	Width:	11.75"	(300 mm)
	Depth:	3.75"	(96 mm)
Construction:	19 gauge, cold-rolled steel		
Transformer:	24 VAC, 30 VA (SS106)		
Storage:	Two 12 VDC batteries (BA001)		
Connection:	Screw terminals (on PK124)		
Mounting:	Surface		



PK124K Battery Backup Kit

Required Components

- CM800 Emergency Call System
- DS100 Door Sentry & Emergency Alarm System
- Tek-ENTRY® Apartment Entry Systems
- or-
- Tek-SAFE® Area of Rescue Assistance System

Optional Components


- Mier BW3000 Offset Cam Lock, or equivalent

www.tektone.com

277 Industrial Park Road • Franklin, NC 28734 • tektone@tektone.net

Phone: (828) 524-9967 • Fax: (828) 524-9968 • Sales: Choose option 2 • Tech Support: Choose option 3
TekTone®'s quality system is registered by UL® to the ISO 9001 standard. (Reference #10001510.)


Figure 1—PK124 to PK151A or PK152 Interconnect Wiring Diagram
Tek-CARE® NC110, NC150 & NC200 Nurse Call Systems; Tek-SAFE® Area of Rescue Assistance System



Notes:

1. Do not connect transformer to PK151A or PK152.
2. Use 24 VAC connections on SS106 transformer.
3. All wiring #18 unless shown.
4. See installation manual for additional PK151A or PK152 connections.

Figure 2—PK124 to PK2019/A or SK109 Interconnect Wiring Diagram
CM800 & CM900 Emergency Call Systems; DS100 Emergency Call System




IL426 PK124 to PK2019 or SK019 Interconnect Diagram REV2 011711 1

Notes:

1. Do not connect transformer to PK2019/A or SK019.
2. Use 24 VAC connections on SS106 transformer.
3. All wiring #18 unless shown.
4. See installation manual for additional PK2019/A or SK019 connections.

Figure 3—PK124 to PK400 Interconnect Wiring Diagram
SM400-series Intercom Systems



IL426 PK124 to PK400 Interconnect Diagram REV1 011711 1

Notes:

1. Do not connect transformer to PK400.
2. Use 24 VAC connections on SS106 transformer.
3. All wiring #18 unless shown.
4. See installation manual for additional PK400 connections.



NC399P Alphanumeric Pager

(CommtechWireless 7900)

IL734 Specification Sheet

Tek-CARE® Nurse Call Systems • Rev. 7 - 11/2012

Features

- Alphanumeric
- 2-line backlit display
- Easy to use large 4-button design
- Configurable for audible or vibration alerts
- Field programmable
- Includes holster belt clip and safety cord

Description

The NC399P alphanumeric pager is a radio pager receiver with a large 2-line display, and is designed for use with TekTone®'s Tek-CARE®500 Wireless Emergency Call System, and TekTone®'s Tek-CARE®400, NC300II, NC200, NC150 and NC110 Nurse Call Systems. The pager stores up to 60 messages. Messages may be selectively deleted or locked into memory. Messages not locked into memory are automatically deleted as new messages arrive (beyond the 60 message limit), with the oldest messages being deleted first. The NC399P may be configured for annunciation via vibration or audible tones, and has one alarm for daily annunciation. It can be programmed to turn on and off at preset times each day.

The NC399P has been pre-programmed to receive call information from TekTone® systems. Its CAP code may be reprogrammed in the field.

Specifications

Dimensions:	Height:	1.9"	(49 mm)
		without holster belt clip	
	Width:	2.7"	(70 mm)
	Depth:	0.7"	(19 mm)
Construction:	Electronic components surface mounted to PCB enclosed in ABS plastic housing		
Operating Freq.:	457.550 Mhz		
Power supply:	1.5 vDC (1 AA cell)		



NC399P Alphanumeric Pager

Applicable Systems

- Tek-CARE®500 Wireless Emergency Call System
- Tek-CARE®400 P5, 400, NC300II, NC200, NC150, NC110 Nurse call systems
- NC365/A/AT Transmitter/Control Unit
- NC366A/AT Transmitter/Control Unit
- NC375 Transmitter Interface

Paging is an ancillary method of call annunciation, and is not intended to replace the primary annunciation provided by the nurse call system.

www.tektone.com

277 Industrial Park Road • Franklin, NC 28734 • tektone@tektone.net

Phone: (828) 524-9967 • Fax: (828) 524-9968 • Sales: Choose option 2 • Tech Support: Choose option 3

TekTone®'s quality system is registered by UL® to the ISO 9001 standard. (Reference #10001510.)



NC377-series Voltage Interface IL680 Specification Sheet

Tek-CARE® NC110, NC150, NC200 & Tek-ALERT® Systems • Rev. 10 - 08/2014

Description

The NC377-series Voltage Interface is used with the Tek-ALERT® Alert Integration Manager to connect switched voltage analog systems to the Tek-ALERT® system. Tek-ALERT® enables such advanced features as consolidated event monitoring, management reporting, and paging for all connected systems. TekTone®'s Tek-CARE® NC110, NC150 and NC200 Nurse Call Systems, as well as a variety of non-TekTone® systems, use the voltage interface to communicate with the Tek-ALERT® system. The voltage interface recognizes three states: off, steady on, and flashing; and can connect systems with input voltage from 5V to 50V, AC or DC.

The NC377-series Voltage Interface can also integrate TekTone®'s Tek-CARE® NC110, NC150 and NC200 Nurse Call Systems directly with the NC369 Paging Transmitter. In this scenario, it connects directly to the nurse call master. The NC377-series Voltage Interface sends nurse call priority information via a serial port to the NC369 Paging Transmitter. The transmitter then transmits this information to the pagers, allowing the staff to distinguish call priority levels.

The voltage interface is available in two configurations: The NC377/32 accommodates 32 inputs, while the NC377/64 handles up to 64 inputs. Connect up to 9 total NC377/64s to the NC369 Paging Transmitter to accommodate up to 576 inputs. (Maximum cable distance from the first NC377-series to the last one is 1,000 feet.)

Specifications

Dimensions: 14"×10"×1.5" (356 mm × 254 mm × 38 mm)
Connections: Screw Terminals: 9-PIN serial connector
Housing: Surface mounted
Voltage: Trigger: 24 volts DC
Operational: 110 VAC Input
Current
Consumption: Less than 450 mA
Messages
Output: 4 per input: 1-Routine, 2-Bath,
3-Emergency 4-Cancel/Reset
BAUD Rate: 9600



NC377-series Voltage Interface
(exterior & interior)

Required Components

NC369 Paging Transmitter
NC399P Pagers
and Tek-CARE® NC110, NC150, NC200 Nurse Call System
or Tek-ALERT® Alert Integration System

Paging is an ancillary method of call annunciation, and is not intended to replace the primary annunciation provided by the nurse call system.

www.tektone.com

277 Industrial Park Road • Franklin, NC 28734 • tektone@tektone.net

Phone: (828) 524-9967 • Fax: (828) 524-9968 • Sales: Choose option 2 • Tech Support: Choose option 3
TekTone®'s quality system is registered by UL® to the ISO 9001 standard. (Reference #10001510.)



NC369 UHF Paging Transmitter

IL670 Specification Sheet

Tek-CARE® Nurse Call Systems • Rev. 10 - 11/2012

Description

The NC369 UHF Paging Transmitter receives alphanumeric messages and patient call information via serial interface and transmits the information directly to pagers. Information is sent to the transmitter at 9600 baud to minimize bottlenecks. Call type and location are determined by the system supplying information to the transmitter. Up to 1000 pager IDs are supported. Battery backup is also available. Programming is simple using a PC or laptop and web based programming screens. Dipswitches are also available on the unit for minimal function onsite setup.

The system's range depends upon the site and the transmitter's location. To determine the exact coverage, conduct a site survey. Estimate a range of up to one mile under ideal conditions. The unit includes a ¼-wave BNC-terminated antenna to provide coverage that meets most requirements. The optional CA133 Antenna is suited for installations that require good building penetration and range. On average, the CA133 Antenna adds about 25% to the transmission footprint, and because it comes with approximately 10 feet of cable, it can be mounted in more locations.

Specifications

Dimensions:	2.7" × 10" × 8" (6.86 cm × 25.4 cm × 20.32 cm)
Construction:	Electronic components mounted on a PCB enclosed in an ABS housing
Connections:	Plug-in RJ45 to DB9
Housing:	Surface or shelf mounted
Frequency:	Programmable for 440–470 MHz
RF Power Output:	Up to 2 watts
Voltage:	12VDC at 2A Input 115 VAC
FCC ID No.:	T5GPT5
IC (Canada):	IC4767A-PT5
BAUD Rate:	9600 serial input, 1200 transmit



NC369 UHF Paging Transmitter

Required Components

- NC399P Pagers
- Tek-CARE®400 P5 Nurse Call System
- or- Tek-CARE®400 Nurse Call System
- or- Tek-CARE®NC300II Nurse Call System
- or- Tek-CARE®NC110/150/200 Nurse Call System (requires NC375 Transmitter Interface)
- or- Tek-CARE®500 Wireless Emergency Call System

Paging is an ancillary method of call annunciation, and is not intended to replace the primary annunciation provided by the nurse call system.

www.tektone.com

277 Industrial Park Road • Franklin, NC 28734 • tektone@tektone.net

Phone: (828) 524-9967 • Fax: (828) 524-9968 • Sales: Choose option 2 • Tech Support: Choose option 3
TekTone®'s quality system is registered by UL® to the ISO 9001 standard. (Reference #10001510.)



SI005 LED Messaging Sign

IL890 Specification Sheet

Tek-CARE®400 P5, Tek-CARE®400 & Tek-CARE®500 Systems • Rev. 3 - 11/2012

Features

- One-line backlit LED display
- Displays 10 to 40 alphanumeric characters per line
- Configurable for audible alerts
- Programmable for up to 6 CAP codes
- Wireless data input
- Receives text messages from either TekTone® system

Description

The SI005 LED Message Sign is a radio paging receiver with a large one-line alphanumeric display. It is designed for use with TekTone®'s Tek-CARE®500 Wireless Emergency Call and Tek-CARE®400 Nurse Call Systems. The SI005 stores up to 24 messages. The messages are deleted after a predetermined length of time. Messages are automatically deleted as new messages arrive (beyond the 24 message limit), with the oldest messages being deleted first. The SI005 may be configured for annunciation via audible tones. The SI005 has been pre-programmed to receive call information from TekTone® systems, and it can be programmed for up to six CAP codes.

Specifications

Dimensions: 28" L x 6.75" H (711 mm x 171 mm)
Weight: 8.4 lbs.
Construction: LED displayboard with RF paging circuitry enclosed in a plastic housing.
Operating Freq.: 457.550 Mhz
Power supply: Input 115 VAC



SI005 LED Messaging Sign

Required Component

NC369 Paging Transmitter

Applicable Systems

Tek-CARE®400 P5 Nurse Call System
Tek-CARE®400 Nurse Call System
Tek-CARE®500 Wireless Emergency Call System

The LED message sign is an ancillary method of call annunciation, and is not intended to replace the primary annunciation provided by the nurse call system.

www.tektone.com

277 Industrial Park Road • Franklin, NC 28734 • tektone@tektone.net

Phone: (828) 524-9967 • Fax: (828) 524-9968 • Sales: Choose option 2 • Tech Support: Choose option 3
TekTone®'s quality system is registered by UL® to the ISO 9001 standard. (Reference #10001510.)